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
1 Fanalysis'e Giriř

Fanalysis Login

Please fill in your credentials to login.

Username

Password

Login

[Lost Password.](#)
Don't have an account? [Sign up now.](#)

Figure 1 Fanalysis Giriř Ekranı

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Fanalysis Sign Up

Please fill this form to create an account.

Username

Password

Confirm Password

The new password must contain at least six characters and categories among the following:

- Uppercase characters (A-Z)
- Lowercase characters (a-z)
- Digits (0-9)
- Special characters (~!@#\$%^&* _+=|(){}[])

Firstname

Lastname

E-Mail

Working Firm

Phone

Already have an account? [Login here.](#)

Figure 2 Fanalysis new user Registration screen

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- Welcome to Fanalysis -

Please select Frontal, Basal or Lateral tab first.

Drag and Drop a file here
 or Double click into area to upload with file selector!
 After upload, Click on an image to use on tab.

Save Frontal image as
Save Frontal
Save All
Load Project
Print

Auto select next option after clicked
 Auto hide lines and points while pointing

Frontal

Lateral

Basal

Operative Plan

Saved Projects

Photographic Analysis

Aesthetic Analysis

Ruler length: mm Supratip (S)

- Ruler Begin (RB)
- Ruler End (RE)
- Vertex (V)
- Trichion (Tr)
- Glabella (G)
- Lateral border of Head - Left (LHL)
- Lateral border of Head - Right (LHR)
- Nasion (N)
- Rhinion (R)
- Endocanthus - Left (ENL)
- Endocanthus - Right (ENR)
- Lateral canthus - Left (LCL)
- Lateral canthus - Right (LCR)
- Widest point of dorsum - Left (WDL)
- Widest point of dorsum - Right (WDR)
- Lower bony vault - Left (XL)
- Lower bony vault - Right (XR)

- Tip (T)
- Tip defining point - Left (DL)
- Tip defining point - Right (DR)
- Alar flare - Left (ALL)
- Alar flare - Right (ALR)
- Subnasale (Sn)
- Labrale superius (LS)
- Stomion (St)
- Menton (Me)
- Medial brow - Left (MBL)
- Medial brow - Right (MBR)
- Lateral brow - Left (LBL)
- Lateral brow - Right (LBR)
- Lateral limbus - Left (LLL)
- Lateral limbus - Right (LLR)
- Medial limbus - Left (MLL)
- Medial limbus - Right (MLR)

Photometric Analysis

Figure 3 Fanalysis logged in screen view

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Figure 4 Fanalysis user update and logout links

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
Figure 5 Fanalysis process menu sections

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
Figure 7 Picture loaded and selected from the loader area

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
Figure 8 Photographic analysis options for Frontal area

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
Figure 9 Working section options selector

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
Figure 10 Frontal/Lateral/Basal Calculations values

Photographic Analysis **Aesthetic Analysis**

Frontal Operative analysis

Face Shape	-		
6:8	<input checked="" type="radio"/> Proportional <input type="radio"/> Disproportional		<input type="checkbox"/> Scar
Horizontal thirds	<input checked="" type="radio"/> Proportional <input type="radio"/> Disproportional		<input type="checkbox"/> Dorsum
Vertical fifths	<input checked="" type="radio"/> Proportional <input type="radio"/> Disproportional		<input type="checkbox"/> Sidewall-left
Facial asymmetry (specify)	<input type="checkbox"/>		<input type="checkbox"/> Sidewall-right
Brow-Tip Aesthetic Line	Right: <input checked="" type="radio"/> Smooth <input type="radio"/> Irregular		<input type="checkbox"/> Radix
	Left: <input checked="" type="radio"/> Smooth <input type="radio"/> Irregular		<input type="checkbox"/> Tip
Skin	Thickness: -		<input type="checkbox"/> Alar-left
	Sebaceous glands: -		<input type="checkbox"/> Alar-right
	Skin discoloration: -		Specify
	Mobility of the skin: <input checked="" type="radio"/> Normal <input type="radio"/> Limited		<input type="checkbox"/> Bony vault Right
	<input type="checkbox"/> Freckles		<input type="checkbox"/> Bony vault Left
			<input type="checkbox"/> Tip Right
			<input type="checkbox"/> Tip Left
			<input type="checkbox"/> Cartilaginous vault Right
			<input type="checkbox"/> Cartilaginous vault Left
			Radix -
			Bony vault -
			Cartilaginous vault -
			Tip -
			Alar axis -
			Alar size -
			Alar symmetry <input type="checkbox"/>
			Alar Flare -
			Infratip lobule -
			Chin width -
			<input type="checkbox"/> Supraalar pinching
			<input type="checkbox"/> Tip bifidity
			<input type="checkbox"/> Parenthesis deformity
			<input type="checkbox"/> Unitip deformity
			<input type="checkbox"/> Bulbosity
			<input type="checkbox"/> Asymmetry
			<input type="checkbox"/> Cephalically oriented lateral crura
			Additional remarks

Figure 11 Aesthetic Analysis section options for Frontal area

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Save Save Operational Plan Save All Load Project Print Auto select next option after clicked Auto hide lines and points while pointing

Frontal Lateral Basal **Operative Plan** Saved Projects

Operative Plan

Projection	Do not disturb
Rotation	Do not disturb 111 mm
Type of procedure	Rhinoplasty
Inferior turbinate surgery	-
Other intranasal surgery	Partial resection
Incisions	Killian
Approach	Delivery
Septum	Cartilage harvest
Dorsum (Profile alignment)	Lateral view: Augmentation Frontal view: Do not disturb
Radix	Do not disturb
Tip	Lateral crura: Do not disturb Reorientation sutures: <input checked="" type="checkbox"/> Tongue-in-Groove <input checked="" type="checkbox"/> Transdomal suture <input checked="" type="checkbox"/> Domal suture Left+Right <input checked="" type="checkbox"/> Septocolumellar suture <input checked="" type="checkbox"/> Interalar suture <input checked="" type="checkbox"/> Intermediate crura sutured together <input checked="" type="checkbox"/> Footplate narrowing suture Other(s): 1) other1 2) other2 3) other3 Medial crura Caudal margin resection
Nasal Packing	Vaseline gauze Splint: Left Merocel: Right Vaseline gauze: Left+Right

Alar Reduction	<input checked="" type="checkbox"/> Nostril <input checked="" type="checkbox"/> Nostril floor <input checked="" type="checkbox"/> Alar lobule																																													
Nasolabial angle	Augmentation																																													
Osteotomies	Medial: Left Lateral: Right Intermediate: Bilateral Transverse: None																																													
Other procedures	<input checked="" type="checkbox"/> Mentoplasty <input checked="" type="checkbox"/> Blepharoplasty <input checked="" type="checkbox"/> Closure of septal perforation Other: other																																													
Miscellaneous	<input checked="" type="checkbox"/> Quilting suture <input checked="" type="checkbox"/> Division of depressor septi muscle <input type="checkbox"/> Vestibular skin resection																																													
Grafts	<table border="1"> <thead> <tr><th>Graft type</th><th>Side</th><th>Material</th></tr> </thead> <tbody> <tr><td><input checked="" type="checkbox"/> Columellar strut</td><td></td><td>Septal Cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Spreader graft</td><td>Left</td><td>Auricular cartila</td></tr> <tr><td><input checked="" type="checkbox"/> Spreader flap</td><td>Right</td><td></td></tr> <tr><td><input checked="" type="checkbox"/> Shield graft</td><td></td><td>Costal cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Tip graft (Peck)</td><td></td><td>Septal Cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Lateral crural graft</td><td>Bilater</td><td>Septal Cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Lateral crural strut graft</td><td>Left</td><td>Septal Cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Alar batten</td><td>Right</td><td>Auricular cartila</td></tr> <tr><td><input checked="" type="checkbox"/> Rim graft</td><td>Bilater</td><td>Costal cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Radix graft</td><td></td><td>Costal cartilage</td></tr> <tr><td><input checked="" type="checkbox"/> Caudal septal extension graft</td><td></td><td>Auricular cartila</td></tr> <tr><td>Other 1)</td><td>other1</td><td></td></tr> <tr><td>Other 2)</td><td>other2</td><td></td></tr> <tr><td>Other 3)</td><td>other3</td><td></td></tr> </tbody> </table>	Graft type	Side	Material	<input checked="" type="checkbox"/> Columellar strut		Septal Cartilage	<input checked="" type="checkbox"/> Spreader graft	Left	Auricular cartila	<input checked="" type="checkbox"/> Spreader flap	Right		<input checked="" type="checkbox"/> Shield graft		Costal cartilage	<input checked="" type="checkbox"/> Tip graft (Peck)		Septal Cartilage	<input checked="" type="checkbox"/> Lateral crural graft	Bilater	Septal Cartilage	<input checked="" type="checkbox"/> Lateral crural strut graft	Left	Septal Cartilage	<input checked="" type="checkbox"/> Alar batten	Right	Auricular cartila	<input checked="" type="checkbox"/> Rim graft	Bilater	Costal cartilage	<input checked="" type="checkbox"/> Radix graft		Costal cartilage	<input checked="" type="checkbox"/> Caudal septal extension graft		Auricular cartila	Other 1)	other1		Other 2)	other2		Other 3)	other3	
Graft type	Side	Material																																												
<input checked="" type="checkbox"/> Columellar strut		Septal Cartilage																																												
<input checked="" type="checkbox"/> Spreader graft	Left	Auricular cartila																																												
<input checked="" type="checkbox"/> Spreader flap	Right																																													
<input checked="" type="checkbox"/> Shield graft		Costal cartilage																																												
<input checked="" type="checkbox"/> Tip graft (Peck)		Septal Cartilage																																												
<input checked="" type="checkbox"/> Lateral crural graft	Bilater	Septal Cartilage																																												
<input checked="" type="checkbox"/> Lateral crural strut graft	Left	Septal Cartilage																																												
<input checked="" type="checkbox"/> Alar batten	Right	Auricular cartila																																												
<input checked="" type="checkbox"/> Rim graft	Bilater	Costal cartilage																																												
<input checked="" type="checkbox"/> Radix graft		Costal cartilage																																												
<input checked="" type="checkbox"/> Caudal septal extension graft		Auricular cartila																																												
Other 1)	other1																																													
Other 2)	other2																																													
Other 3)	other3																																													
Camouflage / Material	<input checked="" type="checkbox"/> Radix: Cartilage chips <input checked="" type="checkbox"/> Middle third: Thinned cartilage <input checked="" type="checkbox"/> Tip: Cartilage dust																																													
Lateral Profile Changes (mm)	Radix: radix Rhinion: rhinion Middle Vault: middle vault Supratip: supratip Tip: tip Subnasale: subnasale																																													
Additional remarks	additional																																													

Figure 12 Fanalysis Operative Planning options

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
Figure 14 Fanalysis saved projects page

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©2022 - Fanalysis.net

Figure 15 Fanalysis Reports example for Frontal section